

# RAMCO Catalogue 2014

Dear Ladies and Gentlemen and  
dear friends of RAMCO.

We are pleased to present our latest RAMCO-Catalogue "March 2014".

Please note, that the articles shown in this catalogue is only part of our  
range. We got more than 10,000 articles numbers in our range.

All deliveries, packaging and clearance from our warehouse is under  
control and supervision by our in-house QualityManagement.

Please direct enquiries, claims and suggestions directly to:


QM department  
ms. Annemette Rasmussen  
Phone: +45 2396 3999  
Telefax: +45 86 401 789  
email: am@ramco.dk

All products is CE-Marked and ALL our suppliers / manufactures is  
ISO 9001:2008 and / or ISO 13485:2003 certified.

No liability for errors or misprints.  
© RAMCO Ltd.  
March 2014

# CONTENTS

|  |         |
|--|---------|
| Diagnostica ( kits, reagents & Point Of Care ) .....         | 3 - 10  |
| Disposables .....  | 11 - 20 |
| Equipment .....  | 21 - 24 |
| Wound Care Products .....                                    | 25 - 28 |
| Instruments .....  | 29 - 32 |
| First-Aid Products .....                                     | 33 - 36 |
| Medical kits and Tailor-Made kits ("Kitting Services") ..... | 37 - 38 |
| Dental products .....  | 39 - 40 |
| Company profile .....  | 41      |



**SECTION: Diagnostica**



### Malaria Test

One Step Malaria P. F. Test is a simple, rapid test for the diagnosis of malaria using serum/plasma samples. One Step malaria P. F. Test detects of malaria parasites and for differentiation of Plasmodium Falciparum (P. F.) from other less virulent malaria species.

One Step Malaria Pf/Pv Test is a lateral flow chromatographic immunoassay for the detection of Malaria histidine-rich protein II (Pf HRP-II) and plasmodium lactate dehydrogenase (pLDH) in vitro assay. It is used for detection of Malaria Plasmodium falciparum and Plasmodium Vivax in human blood. This device is intended to be used as a screening test and as an aid in the diagnosis of infection with Plasmodium.

Registration: CE / WHO

**RAMCO Art.No. 46000**



| Product       | Specimen             | Product format    | Sensitivity                          | Specificity | Reaction time | Shelf life |
|---------------|----------------------|-------------------|--------------------------------------|-------------|---------------|------------|
| Malaria pf    | Serum, Plasma, or WB | Strip or cassette | 20 parasites/ $\mu$ L (WHO Standard) | > 99%       | 15 min        | 24 months  |
| Malaria pf/pv | Serum, Plasma, WB    | Strip or cassette | 20 parasites/ $\mu$ L (WHO Standard) | > 99%       | 15 min        | 24 months  |

### HBV (Hepatitis B) Test

One Step HBV [5 in 1] Test is a rapid, qualitative immunoassay for the determination of hepatitis B markers (HBsAg, HBsAb, HBeAg, HBeAb, HBcAb) in human serum/plasma/whole blood as an aid in the clinical assessment of hepatitis B infection.

Box with 20 test

**RAMCO Art.No. 46005**



### HCV (Hepatitis C) Test

One Step Anti-HCV Test is a rapid direct binding test for the qualitative detection of anti hepatitis C antibodies in human serum/plasma/whole blood as an aid in the clinical assessment of hepatitis C infection. Test results are read visually without any instrument.

Available in Box with 25, 50 or 100 pcs

**RAMCO Art.No. 46008**





### HIV Rapid test

HIV 1+2 Test is a rapid direct binding test for the qualitative detection of antibodies to Human Immunodeficiency Virus Type 1 and Type2 (HIV) in human serum/plasma or whole blood as an aid in the clinical assessment of HIV infection.

This one step test is very sensitive and only takes about 10-15 minutes. Test results are read visually without any instrument.

Test is available in both Strip and Cassette format.

Packed with 25 pcs / box.

**RAMCO Art.No. 46010 (Cassette) og 48012 (Strip)**



### Syphilis Rapid test

Syphilis is a test to qualitatively detect Treponema pallidum antibodies in human whole blood, serum or plasma.

Application: apply sample, result can be read after 10 minutes.

Relative Sensitivity: 99,7 %

Relative Specificity: 99,6 %

Accuracy: 99,7 %

Storage: 2 - 30 °C.

Shelf life: 24 months

**RAMCO Art.No. 46013**



### Strep-A Test

A rapid test for the qualitative detection of group A streptococcal antigens in throat swabs.

The results are easy to read:

2 stripes show a positive result,

1 stripe means negative.

The control stripe indicates that the test has been performed correctly.

Results appear within 5 minutes.

Kit complete - ready to use incl. sample swabs.

Available in both Strip and Cassette format.

Storage: 2 - 30 °C.

**RAMCO Art.No. 48120 (Cassette) og 48121 (Strip)**



### Influenza A+B Test

Influenza Antigen Type A and Type B.

Quick Chromatographic immunoassay to qualitatively detect influenza A- and B-antigenes from nasal swabs.

Result after 10 minutes.

Storage: 2 - 30 °C

Available in 5, 10 or 25 test / box.

**RAMCO Art.No. 46118**





**Tuberculosis test**

This rapid, qualitative chromatographic immunoassay test detects human Anti-TB (Myobacterium tuberculosis) antibodies in serum or plasma.

Conducting the test is simple:

Place the sample onto the test cartridge.

Results in 15 minutes.

Storage at 2 - 30 °C.

Pack available in 3 sizes ( 5, 10 or 25 test).

Contents in box:

Test cartridges, disposable pipettes

**RAMCO Art.No. 46124**



**Gonorrhoeae (NGH) Test**

One Step Neisseria Gonorrhoeae (NGH) Test is a rapid and convenient immunochromatographic in vitro assay. It is used for detection of gonorrhoea antigen in urine and the secretory products from urinary system, as an aid in the diagnosis of gonococcus infection. The test provides a visual, qualitative result, and all positive specimens are advised to be confirmed with other qualified assays.

Sensitivity: >96%

Specificity: >99%

Shelf life: 24 months.

**RAMCO Art.No. 46128**



**Fecal Occult Blood (FOB)**

One Step Fecal Occult Blood (FOB) Test is an immunochemical assay with much higher sensitivity and is designed to be more specific in detecting low levels of human fecal occult blood. It does not require special dietary restriction on patients. One Step Fecal Occult Blood Test is a rapid test that can be performed without the use of an instrument. Test yields a positive result when the concentration of hHb in test specimen is at or above 100ng/ml in Sensitivity and the test has a Specificity >99%. Shelf life: 24 months.

**RAMCO Art.No. 46130**



**hCG Test**

Quick Response hCG One-step chromatographic color-label immunoassay for the qualitative detection of hCG in urine or serum. Results in 5 minutes. Sensitivity 25 mIU/ml. Sample size of 4 drops. Room Temp. storage. Long Shelf Life: 18 months. . Easy to Use: No special training needed. Ability to detect low hCG levels as early as 2 days before expected menses.

Sample = Serum or Urine.

Folio pouch contains one test and one pipette.

Test available as Strip or Cassette format.

Box with 10, 25 or 50 test.

Shelf life: 36 months

**RAMCO Art.No. 48006**






**CRP Test**

Quick Response CRP - C.Reactive Protein.  
 The C-Reactive Protein (CRP) is an acute-phase protein, the concentration of which increases at the presence of bacterial or viral infection. Its detection can give information to help evaluate the type of infection present, for example a streptococcus a infection.  
 Elevated CrP levels in inflammation-free patients indicate an increased risk of infarction.  
 Concentrations up to 10 mg/l are regarded as normal in an adult.  
 Storage: 2 - 30 °c.  
 Box with 20 test.

RAMCO Art.No. 48920



**Chlamydia Test**

Immunological rapid test to qualitatively detect the chlamydia antigen in the female cervical swab and the male urethral swab or urine sample.  
 After 10 minutes a clear colouration denotes a positive or negative result.  
 Cartridge test, sealed singly.  
 Storage: 2 - 30 °c.  
 Sensitivity: >97%  
 Specificity: >99%

RAMCO Art.No. 46133



**HemoCue Products**

Point-of-care medical devices from HemoCue® deliver laboratory accurate results within minutes, enabling faster diagnosis and treatment. HemoCue® solutions are used in numerous healthcare settings, meeting various customer needs, including connectivity, mobility and user-friendliness

| Art.No | Size                           | Pcs/box |
|--------|--------------------------------|---------|
| 40104  | Hemoglobin Cuvettes            | 50      |
| 40109  | Hemoglobin 201 Cuvettes        | 50      |
| 40105  | Glucose Cuvettes               | 50      |
| 40106  | Glucose 201 Cuvettes           | 50      |
| 40103  | HemoCue Cleaner                | 5       |
| 40108  | HemoCue Hb 201 instrument      | 1       |
| 40110  | HemoCue Glucose 201 instrument | 1       |



**Did you find, what you were looking for??**



**Please call us for any inquiries at +45 86 401 789 or e-mail us at ramco@ramco.dk**

**& we will find the products for you at very competitive prices!!!!**



Test Strips Urine

We provide high quality strips for the most important tests, with long expiry dates and a clear colour chart.  
 For economical and immediate diagnosis.  
 Test strips made in Germany by Meditest.  
 Long Shelf Life time: 24 months.

| Art.No | Description | Analyzes  |
|--------|-------------|---|
| 41225  | Combi-2     | Glucose -3  |
| 41213  | Glucose     | Glucose   |
| 41221  | Glucose -3  | Glucose. Ketone. Ascorbic acid  |
| 41215  | Combi-3a    | Glucose. Ascor. Protein. pH   |
| 41233  | Combi-5     | Gluc. Ascor. Blood. Prot. pH.   |
| 41231  | Combi-5N    | Gluc. Ascor. Blood. Prot. pH. Nitrite.  |
| 41237  | Combi-7     | Gluc. Ascor. Ketone. Blood. Nitrite. pH. Protein                                    |
| 41239  | Combi-9     | Gluc. Ascor. Ketone. Blood. Protein. Nitrite. pH. Bilirubin. Urobilinogen           |
| 41240  | Combi-10L   | Gluc. Ascor. Ketone. Blood. Protein. Nitrite. pH. Bilirubin. Urobilinogen Leukocyte |



Test Strips Blood

Test strips for determination of glucose in blood .  
 We stock several brands.

Example - **Ramco art no 41302** Glucose meter set -  
 Alphacheck Professional:

Automatic encoding, result available in just 5 seconds, only a small amount of blood needed (0.5 µl), blood sampling possible at different body locations, Bio-Sensor Technology, 500 memory locations, test strip ejector. Contents: 1 Glucose meter, 10 strips, 10 lancets, 1 lancing device, 2 batteries, etui



**SIEMENS**  
 Siemens Healthcare  
 Diagnostics





**Blood Lancets**

We got several different bloodlancets.  
 Blood lancets for capillary blood collection, Sterile.  
 Advantages: Practically pain-free, ready for use, also usable at the earlobe, single packed sterile, safe.



**Insuline pen & Needle**

Needles and Pens for injecting Insuline.  
 We got several different brands.  
 Fx Pen Needles - Droplet®  
 sterile, compatible with type A insulin pens.



| Art.No | Size                                | Pcs/box |
|--------|-------------------------------------|---------|
| 46522  | BD Micro-Fine Needle 0,33 x 12,7 mm | 100     |
| 46526  | BD Micro-Fine Needle 0,25 x 8 mm    | 100     |
| 46528  | BD Micro-Fine Needle 0,25 x 5 mm    | 100     |
| 46678  | Penfine Needle 31G x 6 mm           | 100     |
| 46679  | Penfine Needle 31G x 8 mm           | 100     |
| 46680  | Penfine Needle 29G x 10 mm          | 100     |



**NovoPen Insulinesystem**



| Article no. | Description                         |
|-------------|-------------------------------------|
| 46521       | NovoFine Needle 28G x 12 mm 100 pcs |
| 46520       | NovoFine Needle 30G x 8 mm 100 pcs  |
| 46525       | NovoFine Needle 31G x 6 mm 100 pcs  |
| 47082       | NovoPen Insulinepen Blue            |
| 47083       | NovoPen Insulinepen Silver          |





**Resistance testing**


**Resistance testing**

Disk (tablet) approximation tests detecting synergy (or antagonism) are useful to screen and detect the mentioned beta-lactamases.

Storage and use at room temperature, means - Ready for instant use, no need for acclimatization to room temperature - No need for special storage or refrigerator facilities - Long shelf life - generally 4 years.

Neo-Sensitabs range of antibacterials includes the following drug class:

- Penicillins
- Beta-Lactam/ Beta-Lactamase inhibitors Cephems
- Cephamycins & Oxacephems
- Penems and Carbapenems
- Monobactams
- Aminoglycosides
- Tetracyclines
- Fosfomicin
- Chloramphenicol & derivates
- Macrolides, Lincosamides & Streptogramins
- Polypeptides
- Sulphonamides and similars
- Nitrofurans
- Lipopeptide
- Ketolide



**Prefilled Formalincontainer**

CELLSTOR POT - 20ml - 10% NEUTRAL BUFFERED FORMALIN is a Prefilled Specimen Container used for safe transport and storage of specimens.

Designed and manufactured to make the lives of your laboratory personnel safer and easier by supplying the CellStor Pot partially filled with buffered formalin. The CellStor pots contain an innovative sealing mechanism, which ensures a reliable leak proof seal. These features ensure that the chance of a specimen pot leaking during transport will be dramatically reduced, so you have peace of mind that a sample will reach its destination safely and correctly fixed.

With it is prefilled with buffered Formalin it also reduces your or your personell exposure to Formalin fumes.

| Art.No | Size  | Pcs/box |
|--------|-------|---------|
| 40360  | 10 ml | 50      |
| 40361  | 20 ml | 25      |
| 40362  | 60 ml | 25      |



## RAMCO profile

**RAMCO** is an international supplier of medical- and laboratory disposables and medical- and laboratory equipment and diagnostics kits.

**RAMCO** deliver medical relief materials for use in emergency situations, as well as for various kinds of humanitarian assistance.

### What **RAMCO** do:

**RAMCO** work with the best international manufactures there manufacture products there conform to internationally Quality Management System (QMS) e.g. CE-marking, ISO 13485, United States QS etc.

**RAMCO** supply products directly to you, priding ourselves on being honest, reliable and flexible.

### The **RAMCO** team:

The **RAMCO** team is a team of well educated professionals and have more than 25 years experience.

**RAMCO** employ both nurses and freelance doctors.

The **RAMCO** team would like to share our knowledge with our clients and assist in effective and efficient solutions for every medical project.

The **RAMCO** team deal with everything from locating your needs to dealing with the paperwork and delivering the goods, where ever they are needed.

### **RAMCO** values:

Being honest, reliable, flexible and understanding are some of our core values.

**RAMCO** are pride to say, that we provide the best possible service to our clients. We are able to respond to our clients needs at a very short notice.

**RAMCO** care about the local people and environment in which our products will be used, as well as the people and environment we live in.

**RAMCO** follow our Quality- and Environmental Policy strictly. We work according to ISO 9001:2008 guidelines.



**Providing Healthcare to Upcomming Economies**



RAMCO HealtCare  
DK-8920 Randers NV  
Denmark  
TLF +45 86 401 789  
E-MAIL [RAMCO@RAMCO.DK](mailto:RAMCO@RAMCO.DK)  
[WWW.RAMCO.DK](http://WWW.RAMCO.DK)